



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Autoclave Steam Indicator Tape

Product Code :
NMALV10037



Description :

Autoclave Steam Indicator Tape Light Yellow indicator tapes has stripes that changes colour from yellow to dark brown/black when passes through steam sterilization Cycle (i.e. 20 minutes of exposure to a temperature of 121 \pm 1/2 C or 5 minutes of exposure to a temperature of 134 \pm 1/2C)

Technical Specification :

Autoclave tape is used as a process indicator in steam sterilizers.
It is very well suited for the clinical and research environments.
It features an aggressive, pressure-sensitive adhesive designed to adhere to a variety of wraps so that the pack stays secure during sterilization and can be used on direct contact with substrate
Tape also features a chemical indicator to verify that it has gone through the sterilization process.



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,