



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Bedside Screen 3 folds

Product Code :
BSIDE-2010002



Description :

Available in

1. Epoxy Coated
2. Stainless Steel

Supplied with spring wire and curtains.

Technical Specification :

Designed in compliance with international quality with cotton standard, our product is widely used in hospitals, hotels, aviation and railways.

The bed side screen consists of 50 mm dia revolving casters to ensure easy mobility.



Medical Hospital Equipments Manufacturers

Website: [SITE](#)