

**Product Name :**  
I.C.U Bed Multifunction Electric Five Function**Product Code :**  
NMHOPB00008**Description :**

I.C.U Bed Multifunction Electric Five Function Multi-Function ICU Bed with electrically operated back rest tilting 0-70°, Foot rest tilting 0-40°, trendelenburg tilting 12°, Reverse trendelenburg tilting 12°, Height Adjustment (710mm to 510mm). Size: L 2200mm x W1050mm x H 510-710mm (Adjustable Height). All functions controlled with LINAK (Denmark) power device (Actuator) In build control panel on both side of Guard rail for Nurse/Patients to adjust bed position easily. Patient Lock-Out system for Nurse Staff at foot board. CPR position button for bed to reach 0° position quickly. Removable High Quality ABS engineering Plastic head & foot panel. Sz 1945mm X 900mm Four easy lifting Guard rail controlled by gas spring. Epoxy coated Mild steel frame work and 4 section perforated top. Built-in emergency battery backup maintain all functions during electricity failure. 150mm dia swivel noiseless castors with Central braking system. Provision for I.V. Rod on all sides of the bed.