



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Labour & Delivery Bed (With Accessories) Size : Body
section: 1000(L) x 800(W) x 750(H)mm

Product Code :
NMIVP20008



Description :

Labour & Delivery Bed (With Accessories) Made in 2 sections, body section with legs mounted on rubber tips, legs section on swivel castors with brakes. Leg section can be completely recessed under the body section. Size : Body section: 1000(L) x 800(W) x 750(H)mm., Leg section: 900(L) x 800(W) x 750(H)mm. Supports for leg holders with solid mechanism to fix leg holders at required height and width. Mattress-pad 2 sections. Supplied with 2 legs holders with canvas straps OR 2 knee crutches (as desired) adjustable in height and width. Epoxy coated tubular steel frame. Mattress-pad density : 21-25 kg/m³. Carrying capacity: 150 - 180 kg. Mattress-pad: Approx. 60(H)mm. Supplied with S.S. bowl. Knockdown construction



Medical Hospital Equipments Manufacturers

Website: SITE