



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Anaesthesia Face Mask Silicon

**Product Code :**  
NMANI10003



**Description :**

Anaesthesia Face Mask Silicon - Suitable use in respiratory care of anaesthesia and resuscitator, Body with antistatic property, Good grip.

Packing: Box of 10 pcs.

Size: All pediatric and adult 0, 1, 2, 3, 4, 5 Anaesthesia Face Mask Silicon, Anaesthesia Equipments Exporters, Anaesthesia Equipments Manufacturers

**Technical Specification :**

The masks incorporate an integral standard connector with a proven anatomical design.

The clear design ensures excellent clarity; this means optimum visibility through the mask.

The smaller sizes have a cylindrical shape to assist in establishing and maintaining an airtight seal with infants and neonates.



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)