



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Safety Box - Cardboard, 5 Litre

**Product Code :**  
NMCBS35104



**Description :**

Safety Box - Cardboard, 5 Litre Volume : 5.0 Litres safety box Capacity : 150 auto disable syringes of 0.5ml with needles Wall thickness : 0.8mm, Size : 305x150x110mm(+/-2mm) Printing : 2 colors standard printing. white /yellow background, printing red and black Diameter of syringe insert hole : 35-40mm



**Medical Hospital Equipments Manufacturers**

Website: SITE,