



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg
& Reverse Trendelenburg 12½°)

Product Code :
BCIVP10001



Description :

Infant Bed / Child Cot / Baby Bassinet (with Trendelenburg & Reverse Trendelenburg 12½°) Size: 870mm L
530mm W 920-1070mm H. Height Adjustable by Gas Spring 920-1070mm. Moulded Plastic Crib
(without joints) with clear visibility for observer. Smooth Rotating Castors 100mm Dia will provide silent
transportation. Brakes on opposite castors. With 12½° Trendelenburg (TR) & Reverse Trendelenburg (RTR)
position



Medical Hospital Equipments Manufacturers

Website: SITE