



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Mobile OT Light

Product Code :
NMOTL10005



Description :

Intensity: 70,000 Lux
Halogen Bulb: 24V - 70W
Reflector: Diachronic
Colour Temp.: 4300K
Temp. Increase at Head Area: 2-4C
Doma Dia: 22in
Finish: Oven Baked Mobile OT Light, Hospital OT Lights & Lamps Exporters, Hospital OT Lights & Lamps Manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE,