

Product Name :
Oxygen Cylinder, Aluminium, Empty

Product Code :
NMANI10049



Description :

Oxygen Cylinder, Aluminium, Empty **SALIENT FEATURES :** Light Weight & Portable, Corrosion Resistant, Non-Magnetic, Easy to handle, Refillable, Suitable for hospitals, Nursing Homes & Home Care

Technical Specification :

Working Pressure: 150 kgf/cm²
Certification: ISO
Gas Capacity: 7 Cubic Meter
Cylinder Material: Aluminum
Usage: Industrial

NAUGRA®

Medical Hospital Equipments Manufacturers

Website: [SITE](#)