

**Product Name :**  
Hospital Plasma Freezer**Product Code :**  
NMBBE10003**Description :**

Plasma freezer includes powerful single stage refrigeration system. Plasma freezer consumes less energy and achieves faster temperature. Microprocessor based controlling system offers temperature accuracy in plasma freezers. Plasma freezer has high and low limit alarming for bag safety. Inbuilt chart recorder in plasma freezer offers better temperature recording. Plasma freezer has innovative coil in shelves design. Hospital Plasma Freezer - Plasma freezer also has rechargeable battery backup for temperature indication. Plasma freezers come in an array of specifications. Stainless steel shelves in plasma freezer facilitate storage. Voltage stabilizer and computer interface features are available in plasma freezer. Hospital Plasma Freezer, Blood Bank Equipments Exporters, Blood Bank Equipments Manufacturers