

Product Name :
Hospital Blood Collection Monitor**Product Code :**
NMBBE10004**Description :**

Hospital Blood Collection Monitors - Blood Collection Monitor is a compact instrument to provide smooth and gentle rocking for homogeneous mixing with anticoagulant without clot formation of blood cells during collection of blood from a donor.

- Volume can be set in 1 ml. increments.
- Provision of pausing collection and of changing programmed volume during pause.
- Micro-controller based program.
- Display of weight and volume.
- Auto tare facility to accurate for the weight of the bag.
- Motor activated clamping at the end of the collection.
- Audio visual alarm to alert in case of any abnormal condition.
- Auto calibration.
- Highly user-friendly system with LCD messages.
- Display of the time taken for collection.
- Diagnostic mode for technical analysis.

Hospital Blood Collection Monitor, Blood Bank Equipments Exporters,
Blood Bank Equipments Manufacturers

Website: SITE,

NAUGRA[®]

Medical Hospital Equipments Manufacturers