



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Suction Tube

Product Code :
NMSUC100P2



Description :

Suction Tube: Silicon or PVC transparent tubing with soft sealing grips. Retain elasticity even after prolonged use. ID =10mm & wall thickness = 2.5 mm. Excellent flexible & heat resistant. Odor less & repeatedly autoclaveable. Transparency of tube allows visibility of flowing fluids Suction Tube, Parts Suction Units Exporters, Parts Suction Units Manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE