

Product Name :
Macintosh Type Curved Laryngoscope Blades - Stainless Steel, POLISHED FINISH

Product Code :
NMLYG10B06



Description :

Macintosh Type Curved Laryngoscope Blades - Stainless Steel, POLISHED FINISH Laryngoscope Blade SS. Blade Size 3. Polished. Autoclavable

Technical Specification :

Macintosh Laryngoscope Blades combine durable stainless steel with large diameter, high performance full integrated fiber optic bundles.

Usage: Hospital

Material: Stainless Steel

Automation Grade: Automatic

Product Type: Laryngoscope

Shape: Curved