

**Product Name :**  
Macintosh Type Disposable Fibre Optic Laryngoscope  
Blade (Plastic heel, SS Spatula)

**Product Code :**  
NMLYG10D02



**Description :**

Macintosh Type Disposable Fibre Optic Laryngoscope Blade (Plastic heel, SS Spatula) Laryngoscope Blade - Plastic heel, SS Spatula

**Technical Specification :**

Laryngoscope set is used for larynx examination, anesthesia, and resuscitation for endotracheal intubation.

Macintosh type curved blades are available in dull finish hence reduce reflections. Plastic handles use type C battery that lasts longer.

Laryngoscopes are simply specialized torches with a handle containing the batteries and a blade which holds the bulb.