



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Laryngoscope Handle with Inbuilt Led Bulb

Product Code :
NMLYG10H07



Description :

Laryngoscope Handle with Inbuilt Led Bulb Autoclavable

Technical Specification :

Light weight construction.
Diamond knurled handle made surface a secure grip in the hand.
Practically wear free, standard holding fixture for the blade.
Solid handle made of chromium plated metal.
Solid metal contact for reliable activation of illumination.
Simple change of batteries at the base of the handle.



Medical Hospital Equipments Manufacturers

Website: [SITE](#)