

Product Name :
Sphygmomanometer Aneroid Round Mobile**Product Code :**
NMDIG10015**Description :**

Sphygmomanometer Aneroid Round Shaped Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Round Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE