



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Wrist And Forearm Splint

Product Code :
NMABW10009



Description :

Wrist And Forearm Splint Size : S, M, L, XL (Left, Right) Feature : Pre shaped splint for immobilization in dorsiflexion. High quality fabric and foam padding for comfort wear. Thumb flap with velcro closure. Application : Strain and sprain in the wrist. Early cast removal. Minor fractures of the wrist and distal forearm



Medical Hospital Equipments Manufacturers

Website: SITE