

Product Name :
Macintosh Laryngoscope set with handle

Product Code :
NMLYG10S02



Description :

Macintosh Laryngoscope set with handle, With 3 S.S. blades. Satin finish all stainless steel.

Custom made combination of sets available.

Macintosh Laryngoscope set with handle,

Technical Specification :

Laryngoscope set is used for larynx examination, anesthesia, and resuscitation for endotracheal intubation.

Macintosh type curved blades are available in dull finish hence reduce reflections. handles use type C battery that lasts longer.

Laryngoscopes are simply specialized torches with a handle containing the batteries and a blade which holds the bulb.

Website: SITE,

NAUGRA[®]

Medical Hospital Equipments Manufacturers