

Product Name :
ULTRASOUND THERAPY UNIT**Product Code :**
PEETE10004**Description :**

ULTRASOUND THERAPY UNIT Solid-state, Dual frequency,Dual Head,Ultrasound Therapy unit.
SPECIFICATIONS: Out-put Frequency : Two,1 MHz & 3 MHz. Treatment Head : Therapy HAMMER HEAD, it combines 2cm square & 5cm square Soundheads in ONE. Out-put Modes : Continous & Pulsed.(Adjustable in 5 % increments). Display : Digital output display. Timer : Digital, Auto-cut & Beep Power : 220v. AC. 50hz

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE,