

Product Name :
TRANSCUTANEOUS NERVE STIMULATOR**Product Code :**
PEETE10007**Description :**

TRANSCUTANEOUS NERVE STIMULATOR Imported, light weight, pocket model, Two channel, Digital Tens, works on three AA battery cells. It utilizes solid state circuit to provide maximum reliability and stability.

SPECIFICATIONS: Six different stimulation modes LCD display & digital operation. Intensity adjustable to 0-9. Two channel output. Preset 10 minutes Timer. Large buttons for easy operation. Maximum output 200 Hz for severe pain

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE