

**Product Name :**  
TILT TABLE**Product Code :**  
PETEE10006**Description :**

TILT TABLE Manually operated Tilt Table comes with Three wide straps to hold the patient. Works to standing position manually with worm gear arrangement. Range of tilt is calibrated from 0-90 deg. The upholstered table top is 61cm wide x 195cm long x 80cm high and is fitted on heavy duty 35mm square steel tube legs. Mounted on four lockable 10cm dia. wheels for mobility. Table comes complete with an Utility Tray & Two Gripping Handles for A.D.L. activities. Frame beautifully finished in oven baked light blue shade

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)