




Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Luer Adapter (MOQ: 100,000 Pcs Assorted Size)

Product Code :
Blood0010



Description :

Luer Adapter (MOQ: 100,000 Pcs Assorted Size) The NET Luer adapter is a sterile ,single use device to be used with safety Blood collection sets or luer needles in routine blood collection .It is manufactured from stainless steel and polypropylene and is fitted with a safety valve at one end.The perforated label not only serves to simplify identification ,but also acts as a seal if integrity



Medical Hospital Equipments Manufacturers

Website: SITE,