

**Product Name :**  
CPAP without Manual Ventilation**Product Code :**  
ICU0004**Description :**

CPAP without Manual Ventilation

**Technical Specification :**

12 Inch large screen with desirable color schemes

Easy to operate

Support adult and pediatric ventilation mode

Multi-mode ventilation function (VCV, SIMV, PSV , PCV, CPAP , Manual SIGH)

Intelligence operating system

Air Compressor Support (optional)

Self-testing and calibration.

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](#)