



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Schiotz Tonometer

**Product Code :**  
OPHTH2010014



**Description :**

Schiotz Tonometer: Highly precise eye tonometer for measuring intra-ocular pressure. Its reliability and durability makes this quality measuring instrument indispensable. High quality agate bearing for an extremely long service life. Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 sub-divisions. Perfect reading of the scale with red pointer. All vital parts are made of stainless steel, other parts are chrome-plated. Complete in a case with velvet-look inserts. Supplied with three weights (5.5 g, 7.5 g, 10 g) and a conversion table. Schiotz Tonometer, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)