

**Product Name :**  
Clear Bottom Round Oil Bath**Product Code :**  
NMCBS11046**Description :**

Clear Bottom Round Oil Bath We offer a wide range of oil baths that are designed to accommodate round bottom flasks with seat cover. Temperature Range: 300 Degree C  
Size: 1 - 10 Ltrs

**Technical Specification :**

These baths are of double walled construction, exterior and interior made of stainless steel 304. We offer these oil baths in various specifications thus meeting the varied needs of different industries.

**Temperature Range: 300 Degree C****Size: 1 - 10 Ltrs**

- A clear bottom. A uniquely designed Bath. Ensures no heater burn out due to lack of water.
- An absolute over-temperature alarm with auto cut off.
- Double walled, internal and external chamber full stainless steel construction.
- A floating overshoot alarm set 5°C above set point.
- Absolute over temperature alarm - audio - with auto cut off.
- End cycle visual indication with heater cut off.

**Clear Bottom Round Oil Bath Exporters, Clear Bottom Round Oil Bath Manufacturers India China**

---

## Medical Hospital Equipments Manufacturers

Website: [SITE](#).

**NAUGRA<sup>®</sup>**