

Product Name :
Semi Automatic Rotary Microtome**Product Code :**
NMCBS13010**Description :**

Semi Automatic Rotary Microtome We are offering Semi Automatic Rotary Microtome catering to the needs and requirements of our esteemed customers. This equipment is featured with a standard drive system for smooth sectioning and trimming. We have also installed LCD screen displays section, various trimming thickness range and number of sections in this microtome. There is a lockable hand wheel, interchangeable specimen clamp along with a safety alarm system and huge volume waste tray for easy operation. Semi Automatic Rotary Microtome calibrated controls assure an accurate positioning to make re-cuts possible. Specifications: · Section thickness setting range: 0.25-60 μm · Trimming thickness setting Range: 0.25-60 μm · Setting values: 0.25-1 μm increment step 0.25 μm · 1-10 μm increment step 1 μm · 10-20 μm increment step 2 μm · 20-60 μm increment step 5 μm · Minimum setting value: 0.25 μm Semi Automatic Rotary Microtome, Microtome Instrument Exporters, Microtome Instrument Manufacturers