

Product Name :
Stainless Steel Ultrasonic Cleaner

Product Code :
NMCBS17001



Description :

Stainless Steel Ultrasonic Cleaner Tender Suppliers Africa and Stainless Steel Ultrasonic Cleaner Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Stainless Steel Ultrasonic Cleaner Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Stainless Steel Ultrasonic Cleaner China, Stainless Steel Ultrasonic Cleaner Kenya.

Stainless Steel Ultrasonic Cleaner Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Stainless Steel Ultrasonic Cleaner Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the

Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point), Suriname (Paramaribo)

Technical Specification :

These Ultrasonic cleaners are used in most of the Indian hospitals and Dispensaries.

Specifications Ultrasonic Cleaner Tank Dim.:20cm x 15cm x 7cm liter, Cap.: 1.8 Ltrs U/S

Power: 50 Watts U/S, Freq.: 33 khz Input:180, Vac: 270 Vac

S.M.P.S. ensures smooth working of cleaners within the voltage variation of 180-270 vac

Stainless Steel Ultrasonic Cleaner, Ultrasonic Cleaner Exporters, Ultrasonic Cleaner Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)