



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Sphygmomanometer Aneroid Round Shaped - Mobile

Product Code :
NMCBS22002



Description :

Sphygmomanometer Aneroid Round Shaped - Mobile ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Round Shaped - Mobile, OPD Instrument Exporters, OPD Instrument Manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE