

**Product Name :**  
Suction Apparatus Trolley Model**Product Code :**  
NMCBS24003**Description :**

Suction Apparatus Trolley Model Ideal & suitable for Hospitals, Nursing Homes and Clinics · Cabinet made of mild steel sheet with epoxy powder coated · ¼ HP heavy duty motor. Oil immersed vacuum pump · Two wide mouth glass bottles of approx. 1600 ml. each · Synthetic air tight rubber lids. Overflow cut-off device · Non-collapsible PVC Tubing · Vacuum-700 +/- 10mm Hg · Trolley mounted on friction free wheels Suction Apparatus Trolley Model, Suction Apparatus Exporters, Suction Apparatus Manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)