



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Horizontal Cylindrical Sterilizer

Product Code :
NMCBS25002



Description :

Horizontal Cylindrical Sterilizer Steam Working Pressure 1.26 kgf/cm² (2.2 kgf/cm² in case High speed sterilizer
Steam Working Temperature 121oC Operating Voltage 240, single phase, AC supply, 50 Hz · Sterilization
Period: 45 to 50 minutes · Radiant locking system, cold water drainage system · Steam Exhaust: 5 to 7 provided
with timer from 1 to 60 Min with double safety valve, steam releasing valve, pressure gauge for the pressure
inside the chamber, and one pressure gauge to show the pressure of outer chamber both the chamber made of
stainless steel lid of the unit made of heavy duty gun metal piece for extra safety, gasket, water level indicator to
see the water inside, heavy duty power plug with socket provided with the unit, accessories: one spare element,
and two gasket , our common size 400x1200 mm Horizontal Cylindrical Sterilizer, Medical Sterilizers Exporters,
Medical Sterilizers Manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)