

Product Name :
Multichannel Patient Monitoring System**Product Code :**
NMCBS28002**Description :**

Multichannel Patient Monitoring System We are engaged in manufacturing, supplying and exporting a wide range of Multichannel Patient Monitoring System. All our products are designed in compliance with international quality norms & standards, using premium raw material and latest technology. In addition to this, we check all the products on various parameters to assure flawlessness. Other information: · 10.4" Color TFT display, with maximum 4 waveforms · Parameters including ECG/RESP, SpO2, NIBP, TEMP and pulse rate · Maximum 72-hour graphic and tabular trends of all parameters · Suitable for adult, pediatric and neonatal patient Multichannel Patient Monitoring System, medical monitoring system exporters, medical monitoring system equipments manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE