



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Dressing Drums - Stainless Steel

Product Code :
NMALV10006



Description :

Dressing Drums - Stainless Steel

Technical Specification :

High-quality material is used for engineering

Highly reliable and durable

Can be used to store sterile bathrobes, dressings, and operational linen etc.

VMS dressing drums with self-locking systems in order to give extra strength for long lasting performance.

Our products are highly durable and corrosion resistance that is ensured by heavy stainless steel material.



Medical Hospital Equipments Manufacturers

Website: [SITE](#)