

Product Name :
Drager Savina**Product Code :**
NMCBS42011**Description :**

Drager Savina Ventilation modes - IPPV (CMV), IPPV Assist (CMV Assist) · SIMV, SIMVASB (SIMV/PS) · CPAP, CPAPASB (CPAP/PS) · BIPAP1, (PCV+) (optional), BIPAP1)ASB (PCV+/PS) (optional) · Enhancements - NIV – Non Invasive Ventilation with optimized alarm system and automatic leakage compensation (optional) · Auto Flow®- Automatic adaptation of the inspiratory flow in volume orientated ventilation modes (optional) · LPO - Low Pressure Oxygen. Independent oxygen supply, e.g. with an O2 concentrator (optional) · Graphic screen - Advanced ventilation monitoring (optional) · Nurse call - Connection for transmitting alarm signals to a central Drager Savina, Medical Ventilator exporters, Medical Ventilator manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE