



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
X-Ray Machine Mobile Portable 100 MA

Product Code :
NMCBS47009



Description :

X-Ray Machine Mobile Portable 100 MA X-Ray Machine Mobile Portable 100 Ma. have Semi Feather Touch Control Panel X-Ray Machine Mobile Portable 100 Ma. have Electronic Solid State Timer With 24 Steps X-Ray Machine Mobile Portable 100 Ma. have Shortest exposure time 10ms X-Ray Machine Mobile Portable 100 Ma. have Independent KVp, mA & mA Displays X-Ray Machine Mobile Portable 100 Ma. have Line Voltage/mA Meter for Voltage Adjustment X-Ray Machine Mobile Portable 100 Ma. have X-Ray Exposure on Audible tone X-Ray Machine Mobile Portable 100 Ma. have Mains ON/OFF Feather Touch Membrane Button X-Ray Machine Mobile Portable 100 Ma. have Overload Function To Protect the tube X-Ray Machine Mobile Portable 100 Ma. have Counter Balanced Mechanism & 3 Wheel design for max. Mobility & manoeuvr ability X-Ray Machine Mobile Portable 100 Ma. have Horizontal and Vertical Travels of arms For Better Positioning X-Ray Machine Mobile Portable 100 Ma. have Collimator Operation from Control Panel



Medical Hospital Equipments Manufacturers

Website: SITE