

**Product Name :**  
FTIR Spectrometers

**Product Code :**  
NMCBS54001



#### Description :

FTIR Spectrometers Tender Suppliers Africa and FTIR Spectrometers Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, FTIR Spectrometers Manufacturers, and Hospital Laboratory Instruments Suppliers from India, FTIR Spectrometers China, FTIR Spectrometers Kenya.

FTIR Spectrometers Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), FTIR Spectrometers Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),

---

Suriname (Paramaribo)

**Technical Specification :**

FTIR Spectrometer · New type cube-corner Michelson interferometer features smaller size and more compact structure, providing higher stability and less sensitive to vibrations and thermal variations than conventional Michelson interferometer · Fully sealed damp and dust proof interferometer, adopting high performance, long lifetime sealing material and desiccator, ensures higher adaptability to the environment and increases accuracy and reliability in operation. Viewable window for silica gel enables easy observation and replacement · Isolated IR source and large space heat dissipation chamber design provides higher thermal stability. Stable interference is obtained without the need of dynamic adjustment · High intensity IR source adopts a reflex sphere to obtain even and stable IR radiation · Cooling fan stretch suspending design ensures good mechanical stability · Super wide sample compartment provides more flexibility to accommodate various accessories · The application of programmable gain amplifier, high accuracy A/D converter and embedded computer improves the accuracy and reliability of the whole system · The spectrometer connects to PC via a USB port for automatic control and data communication, fully realizing plug-and-play operation · Compatible PC control with user friendly, rich function software enables easy, convenient and flexible operation. Spectrum collection, spectrum conversion, spectrum processing, spectrum analyzing, and spectrum output function etc. can be performed · Various special IR libraries are available for routine search. Users can also add and maintain the libraries or set up new libraries by themselves · Accessories such as Defused/Specular Reflection, ATR, Liquid cell, Gas cell, and IR microscope etc can be mounted in the sample compartment Specifications: Spectral Range 7800 to 350 cm<sup>-1</sup> Resolution Better than 0.85cm<sup>-1</sup> Wave Number ±0.01cm<sup>-1</sup> Scanning Speed 5-step adjustable for different applications Signal to noise ratio better than 15,000:1 (RMS value, at 2100cm<sup>-1</sup>, resolution: 4cm<sup>-1</sup>, detector: DTGS, 1 minute data collection) Beam splitter Ge coated KBr Infrared Source Air-cooled, high efficiency Reflex Detector DTGS Data system Compatible computer Software FT-IR software contains all routines needed for basic spectrometer operations, including library search, quantitation and spectrum export IR Library 11 IR libraries include FTIR Spectrometer, FTIR Spectrometer exporters, FTIR Spectrometer manufacturers

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)