



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Luxury Electric Blanket

**Product Code :**  
NMCBS93003



**Description :**

Luxury Electric Blanket The blanket is provided of thermal micro chip cutouts which helps in uniform distribution of heat in whole blanket · In luxury range, the size of single bed is of king size i.e. 36 x 60 inches approx. and double beds is of king size i.e. size is 72 x 60 inches resulting in better coverage of bed · The outer shell material used is of high quality which give cozy feel while you sleep on it Technical Specifications: Features Single bed Double bed King Size (Approx.): 36 x 60 inches. 72 x 60 inches Wattage (Approx.): 80 watts. 150 watts. Voltage: (Operable At Low Voltages Too) 220 - 250 volts (AC) 220 - 250 volts (AC) Frequency: 50-60 Hz 50-60 Hz Luxury Electric Blanket, Heating Pads exporters, Heating Pads manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)