



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Syringe Pump Stand

Product Code :
NMCBS109019



Description :

Syringe Pump Stand * Syringe pump stand is with 4 adjustable shelves. * These shelves helps to hold syringe pumps. * With IV stand on top. * Five 50mm double wheel castors with brake. * Electrical connection is provided at the bottom of the stand. * The Complete stand is mobile on 50mm dia twin castors Syringe Pump Stand, Hospital Trolley exporters, Hospital Trolley manufacturers



Medical Hospital Equipments Manufacturers

Website: [SITE](#)