



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Standard Hospital Bedside Locker

**Product Code :**  
NMCBS110002



**Description :**

Standard Hospital Bedside Locker ÷½ Overall size : 40L x 40W x 82H cms ÷½ Machine pressed sheet box with door of heavy guage ÷½ Assembled to CRC tubular legs of 2.5 cms. dia. fitted with PVC stumps ÷½ Stainless steel polished top of 22G ÷½ Pre-treated and epoxy powder coated Standard Hospital Bedside Locker, Hospital Lockers exporters, Hospital Lockers manufacturers



**Medical Hospital Equipments Manufacturers**

Website: SITE,