



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
BED SIDE LOCKERS

Product Code :
NMCBS110007



Description :

BED SIDE LOCKERS Bed side lockers made of high quality sheet metal. 1 1/2 SS top. 1 1/2 Hospital bed side lockers Provided with one storage cabinet. 1 1/2 Bed side lockers optionally available with single draw under the top and space for keeping utilities. 1 1/2 Standard dimension: 160mm L x 160mm W x 320 mm H. BED SIDE LOCKERS, Hospital Lockers exporters, Hospital Lockers manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE