

**Product Name :**  
Plastic Sprayed Bed Side Lockers**Product Code :**  
NMCBS110009**Description :**

Plastic Sprayed Bed Side Lockers The Hospital Plastic Sprayed Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Plastic Sprayed Bedside Locker has a stainless steel top and mounted on two swivel castors. Hospital ABS Bedside Locker. The Hospital ABS Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Bedside Locker has a stainless steel top and mounted on two swivel castors. Plastic Sprayed Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers