



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Stainless Steel Surface Bed Side Lockers

Product Code :
NMCBS110010



Description :

Stainless Steel Surface Bed Side Lockers The Hospital Stainless Steel Surface Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Stainless Steel Surface Bedside Locker has a stainless steel top and mounted on two swivel castors. Stainless Steel Surface Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE