



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Lamp with 4 Reflector

Product Code :
NMCBS113002



Description :

Ceiling Lamp with 4 Reflector & Dome Head: 515mm Dia. & Halogen Bulbs :4nos. 12v x 50 watt & Light intensity at 1 mt.: 70,000 Lux Intensity Control: 4 Stage & Height Adjustment: 550 mm & Action Radius: 1550mm & Possible Movements: Radial Angular & Axial Color & Temperature: 4200 K & 300K & Temp.rise in field: 3-6 c from Amb.Temp & Low voltage unit with CVT: 1 No. (C.V.T optional) & Input Supply Voltage: 220 Volt A.C & 10% 50Hz & 2% Ceiling Lamp with 4 Reflector, Operation Theater Lights exporters, Operation Theater Lights manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE