

Product Name :
Operation Lights**Product Code :**
NMCBS113005**Description :**

Operation Lights A) Dome Head :715mm Dia. · Halogen Bulbs: 7nos. 12v x 50 watt · Light intensity at 1 mt.: 1,20,000 · Lux Intensity Control: 4 Stage · Height Adjustment: 500mm Action · Radius: 1850mm · Possible Movements: Radial, Angular & Axial Color Temperature :4200 K $\pm$ 300K · Temp.rise in field: 4° - 8° c from Amb.Temp · Low voltage unit with CVT: 1 No. (C.V.T optional) · Input Supply Voltage: 220volt A.C $\pm$ 10% 50Hz $\pm$ 2% B) Dome Head :515mm Dia. · Halogen Bulbs: 4nos. 12v x 50 watt · Light intensity at 1 mt.: 70,000 · Lux Intensity Control: 4 Stage · Height Adjustment: 550mm Action · Radius: 1550mm · Possible Movements: Radial, Angular & Axial Color Temperature :4200 K $\pm$ 300K Temp.rise in field :3°-6° c from Amb.Temp · Low voltage unit with CVT:1 No. (C.V.T optional) · Input Supply Voltage: 220volt A.C $\pm$ 10% Operation Lights, Operation Theater Lights exporters, Operation Theater Lights manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE