

**Product Name :**  
Dental Surgical Microscope

**Product Code :**  
NMCBS131002



#### Description :

Dental Surgical Microscope Owing to our experience, we have been successful in catering to the requirements of clients by offering quality Dental Surgical Microscope. All the products are widely acclaimed by the customers for their excellent performance, easy installation, low maintenance and longer service life. Apart from this, we provide customization as per the specifications laid down by the customers. Specification: Binocular Tube: 45° 1/2, F=170mm Magnification: 0.4x, 0.63x, 1.0x, 1.6x & 2.5x with five step magnifier Eye Piece: WF 10x/16mm with eye guards Working Distance: F=200mm objective lens Focusing: Motorized fine focusing Light Source: 15v, 150w twin halogen reflector lamp Built in Filters: Green & Cobalt Blue Vertical Movement of Arm: 600mm Microscope Carriers: XYZ Coupling or Z-axis Ophthalmic carrier Optional: Assisto scope, Xpert Vu Camera Module, Double Beam Splitter, 0-210° 1/2 tiltable ergo tube, eye piece WF 12.5x/16mm Dental Surgical Microscope, Dental Microscopes exporters, Dental Microscopes manufacturers

#### Technical Specification :

Pure Cut gives the surgeon powerful and smooth cutting action.  
Coagulation is suitable for contact and non-contact operation.  
Bipolar output is suitable for gentle bipolar effect in all types of surgeries.  
Separate Monopolar and Bipolar R.F. output system allow two surgeons to work separately.  
Absolutely Safe as Isolated RF Output provides direct current return to the generator significantly reducing the risk of alternate site burns.  
Power-Limiting Output with Automatic Compensation Circuit provides more consistent performance over wide range of tissue impedance.



## Medical Hospital Equipments Manufacturers

**Website:**Medical Hospital Equipments Manufacturers,