



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Diagnostic Sets

Product Code :
NMCBS137013



Description :

Diagnostic Sets Diagnostic Set Includes: - One Ophthalmoscope Head with 6 different apertures micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter. 28 focusing lenses with a range of -25 to +40 diopters. Rubber brow rest prevents scratching of eyeglasses. Otoscope Head with 3x magnification lens fitted with 3.5v Halogen Lamp. One Halogen Otoscope Head with Reusable Ear Specula 3mm & 4mm, One Halogen Ophthalmoscope Head with 6 Apertures, One Spreadable Nasal Speculum, One Bent Arm Illuminator, 2 Laryngeal Mirrors 4 No. & 5 No., One Tongue Depressor, One Tongue Depressor Holder, One 9V cell Dry Battery Handle. Set complete in a Case. Diagnostic Sets, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,