



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Indirect Ophthalmoscope

Product Code :
NMCBS137017



Description :

Indirect Ophthalmoscope Halogen Binocular Indirect Ophthalmoscope is an excellent instrument for all standard applications, as well as with uncooperative pediatric patients & in cases where the papillary diameter is less than 4mm. The need for pupillary dilatation is virtually obviated when using the SPF-capability of this dual purpose clinical instrument. Indirect Ophthalmoscope, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,