



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)  
Phone: 9896600003

**Product Name :**  
Keratometer

**Product Code :**  
NMCBS137020



**Description :**

Keratometer Keratometer is the world wide standard for manual measurement of corneal curvatures. Measures concave surfaces for contact lens fitting. Both major and minor meridians can be measured once the major axis is located. Easy elevation and lowering allows quick patient positioning. Other features include engraved focusing scale, coincidence focusing system, and single position measurement in both meridians. Keratometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Medical Hospital Equipments Manufacturers**

Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)