



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Lensometer

**Product Code :**  
NMCBS137022



**Description :**

Lensometer Manual Lensometer combines the traditional features you would expect, with innovative advantages, such as LED illumination and cord-free battery power. The cool-running, energy-efficient LED light source provides a clear green illumination target, and virtually maintenance-free operation. Battery operated and completely cordless, the KL-18.010 Manual Lensometer gives you the freedom of mobility without any clutter. The KL-18.010 low-power consumption and auto sleep feature ensure long battery life. Lensometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,