



Email : [sales@naugramedical.com](mailto:sales@naugramedical.com)

Phone: 9896600003

**Product Name :**  
Schiotz Tonometer

**Product Code :**  
NMCBS137026



**Description :**

Schiotz Tonometer Tonometer according to Prof. Schiotz for measuring the Intraocular pressure. All main parts are made of High Quality Brass and are chrome plated. Red needle enables easy reading of the scale. The instrument comes complete with 3 weights (5.5 gram, 7.5gram, 10 gram) and conversion table in Blue Presentation case. Schiotz Tonometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,