



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Schiotz Tonometer

Product Code :
NMCBS137026



Description :

Schiotz Tonometer Tonometer according to Prof. Schiotz for measuring the Intraocular pressure. All main parts are made of High Quality Brass and are chrome plated. Red needle enables easy reading of the scale. The instrument comes complete with 3 weights (5.5 gram, 7.5gram, 10 gram) and conversion table in Blue Presentation case. Schiotz Tonometer, Ophthalmic Equipment, Ophthalmic Equipment manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE